



Certificate of Analysis

Product Name: TG
Product: Purified Recombinant Protein
Lot Number:
Manufactured:

Catalog No.: TP701255
Concentration :
EXP date:

Protein Purification

Starting Material: Culture supernatant
Buffer/Components: PBS.
Storage Conditions: -80°C

Test:	Method:	Document:	Results:
Protein Concentration	Bradford protein assay method		
Purity	SDS-Page (reduced)	G3302	

Certified by:

Date:

No.168 Mashang Meiliang Road | Wuxi , Jiangsu province, China | 0510-81830666

www.origene.com

custsupport@origene.com