

**Product Name:** IL6 Mouse monoclonal antibody, clone 84B7  
**Product:** Purified Monoclonal Antibody **Catalog No.:** D624  
**Lot Number:** **Concentration:**  
**Production**  
**Expiring Date:**

**Antibody Purification**

**Starting Material:** cell culture supernatant  
**Axial Flow Column:** Protein A Sepharose  
**Buffer/Components:** PBS (pH7.3), 0.1%NaN3  
**Storage Conditions:** Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	A280	G3476	
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOSTrip™	G3346	
Application	Elisa	Datasheet	

**Prepared By:**

Quality Control-Date

**Reviewed By:**

Quality Assurance-Date

No.168 Mashan Melliang Road | Wuxi, Jiangsu Province, China |0510-81830666

[www.origene.com](http://www.origene.com) [custsupport@origene.com](mailto:custsupport@origene.com)