

**Product Name:** NGAL mouse monoclonal antibody, clone 10F11  
**Product:** Purified Monoclonal Antibody **Catalog No.:** DM0031-B2  
**Lot Number:** **Concentration:**  
**EXP Date:**

**Antibody Purification**

**Starting Material:** cell culture supernatant  
**Axial Flow Column:** Protein G Sepharose  
**Buffer/Components:** PBS (pH7.3), 0.1%NaN3  
**Storage Conditions:** Store at 2-8°C ,Do Not freeze.

<b>Test:</b>	<b>Method:</b>	<b>Document:</b>	<b>Results:</b>
Protein Concentration	BCA	G3576	
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOStrip™	G3346	
Application	Elisa	Datasheet	

**Prepared By:****Reviewed By:**

No.168 Mashan Meiliang Road | Wuxi, Jiangsu Province, China |0510-81830666

[www.origene.com](http://www.origene.com)[custsupport@origene.com](mailto:custsupport@origene.com)