

Product Name: Sheep Anti-human Vitamin D monoclonal antibody ,clone 8F2
Product: Purified Monoclonal Antibody **Catalog No.:** DM0047-B1
Lot Number: **Concentration:**
EXP Date:

Antibody Purification

Starting Material: cell culture supernatant
Axial Flow Column: Protein G Sepharose
Buffer/Components: PBS (pH7.3), 0.1%NaN3
Storage Conditions: Store at 2-8°C ,Do Not freeze.

| Test: | Method: | Document: | Results: |
|-----------------------|--------------------|-----------|----------|
| Protein Concentration | BCA | G3576 | |
| Purity | SDS-Page (reduced) | G3302 | |
| Isotype | ISOSTrip™ | G3346 | |
| Application | ELISA | Datasheet | |

Prepared By: _____
Quality Control-Date

Reviewed By: _____
Quality Assurance-Date