

Product Name: Mouse Anti-human sflt1 monoclonal antibody , clone 4F5
Product: Purified Monoclonal Antibody **Catalog No.:** DM0058-B1
Lot Number: V0101 **Concentration:** 5.35mg/ml
EXP Date: June.2027

Antibody Purification

Starting Material: cell culture supernatant
Axial Flow Column: Protein G Sepharose
Buffer/Components: PBS (pH7.3), 0.1%NaN3
Storage Conditions: Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	BCA	G3576	5.35 mg/mL
Purity	SDS-Page (reduced)	G3302	> 95%
Isotype	ISOSTrip™	G3346	Mouse IgG1
Application	ELISA	Datasheet	PASS

Prepared By: _____
Quality Control-Date

Reviewed By: _____
Quality Assurance-Date