



Certificate of Analysis

Product Name: TP Protein
Product: Purified Recombinant Protein
Lot Number:
Manufactured:

Catalog No.: DP0018-P1
Concentration :
EXP date:

Antibody Purification

Starting Material: cell culture supernatant
Buffer/Components: 20mM PB, pH7.4, 150mM NaCl, 3mM EDTA
Storage Conditions: -20°C, Avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Protein Concentration	BCA methods		
Purity	SDS-Page (reduced)	G3302	

Prepared By:
Quality Control-Date

Reviewed By:
Quality Assurance-Date

No.168 Mashang Meiliang Road | Wuxi , Jiangsu province, China | 0510-81830666

www.origene.com

custsupport@origene.com