



## Certificate of Analysis

**Product Name:** TP Protein  
**Product:** Purified Recombinant Protein  
**Lot Number:**  
**Manufactured:**

**Catalog No.:** DP0018-P2  
**Concentration:**  
**EXP date:**

### Antibody Purification

**Starting Material:** cell culture supernatant  
**Buffer/Components:** 20mM PB, pH7.4, 150mM NaCl, 3mM EDTA  
**Storage Conditions:** -20°C, Avoid repeated freeze-thaws.

| Test:                 | Method:            | Document: | Results: |
|-----------------------|--------------------|-----------|----------|
| Protein Concentration | BCA methods        |           |          |
| Purity                | SDS-Page (reduced) | G3302     |          |

**Prepared By:**  
Quality Control-Date

**Reviewed By:**  
Quality Assurance-Date

No.168 Mashang Meiliang Road | Wuxi , Jiangsu province, China | 0510-81830666

[www.origene.com](http://www.origene.com)

[custsupport@origene.com](mailto:custsupport@origene.com)