



Certificate of Analysis

Product Name: Strep A protein
Catalog No.: DP0038-P1
Lot Number:
Manufactured:
EXP date:

Concentration:

Antibody Purification

Buffer/Components: PBS (pH7.3), 0.1%NaN3
Storage Conditions: Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Application	ELISA	Datasheet	

Prepared By:

Quality Control-Date

Reviewed By:

Quality Assurance-Date

No.168 Mashan Meiliang Road | Wuxi, Jiangsu Province, China |0510-81830666

www.origene.com

custsupport@origene.com