

Product Name: Human IL5
Product: Purified Recombinant Protein
Lot Number:
Manufactured:

Catalog No.: TP750229
Concentration :
EXP date:

Protein Purification

Starting Material: Lysate
Buffer/Components: 50mM Tris, pH8.0, 500 mM NaCl, 10% glycerol, 1% Sarkosyl.
Storage Conditions: -80°C

Test:	Method:	Document:	Results:
Protein Concentration	Microplate BCA method		
Purity	SDS-Page (reduced)	G3302	

Certified by: _____ **Date:** _____