



Certificate of Analysis

PRODUCT CODE: D1930M000-AM

DESCRIPTION: TNT Monoclonal A1

LOT NUMBER:

PRODUCT FORM: Liquid in 20 mM Sodium phosphate, 150 mM Sodium chloride, pH 7.2

BIOLOGICAL SOURCE: Mouse

SPECIFICITY: Reacts with TNT

FILTRATION: 0.22 μ m

PURIFICATION PROCESS: Hydrophobic Interaction Chromatography

PROTEIN CONCENTRATION:

PRESERVATIVES: 0.05% NaN₃

STORAGE: 4°C

pI RANGE: 7.2 – 7.5

Quality Assurance

Date