

Product Name: Human Anti-human SARS-CoV-2 N monoclonal antibody ,clone 7B9
Product: Purified Monoclonal Antibody **Catalog No.:** D632
Lot Number: **Concentration:**
EXP Date:

Antibody Purification

Starting Material: cell culture supernatant
Axial Flow Column: Protein G Sepharose
Buffer/Components: PBS, pH 7.2
Storage Conditions: Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	BCA	G3576	
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOSTrip™	G3346	
Application	ELISA	Datasheet	

Prepared By: _____
Quality Control-Date

Reviewed By: _____
Quality Assurance-Date