

**Product Name:** Mouse Anti-human NP(Influenza A) monoclonal antibody ,clone 5A9  
**Product:** Purified Monoclonal Antibody **Catalog No.:** D675  
**Lot Number:** V0101 **Concentration:** 5.98mg/ml  
**EXP Date:** June.2027

### Antibody Purification

**Starting Material:** cell culture supernatant  
**Axial Flow Column:** Protein G Sepharose  
**Buffer/Components:** PBS (pH7.3), 0.1%NaN3  
**Storage Conditions:** Store at 2-8°C ,Do Not freeze.

| Test:                 | Method:            | Document: | Results:   |
|-----------------------|--------------------|-----------|------------|
| Protein Concentration | BCA                | G3576     | 5.98 mg/mL |
| Purity                | SDS-Page (reduced) | G3302     | > 95%      |
| Isotype               | ISOSTrip™          | G3346     | Mouse IgG1 |
| Application           | ELISA              | Datasheet | PASS       |

Prepared By: \_\_\_\_\_  
Quality Control-Date

Reviewed By: \_\_\_\_\_  
Quality Assurance-Date