

**Product Name:** Mouse Anti-human NP(Influenza A) monoclonal antibody ,clone 8D10  
**Product:** Purified Monoclonal Antibody **Catalog No.:** D676  
**Lot Number:** V0101 **Concentration:** 5.62mg/ml  
**EXP Date:** June.2027

### Antibody Purification

**Starting Material:** cell culture supernatant  
**Axial Flow Column:** Protein G Sepharose  
**Buffer/Components:** PBS (pH7.3), 0.1%NaN3  
**Storage Conditions:** Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	BCA	G3576	5.62 mg/mL
Purity	SDS-Page (reduced)	G3302	> 95%
Isotype	ISOSTrip™	G3346	Mouse IgG1
Application	ELISA	Datasheet	PASS

Prepared By: \_\_\_\_\_  
Quality Control-Date

Reviewed By: \_\_\_\_\_  
Quality Assurance-Date