

**Product Name:** HIV-1 gp41 protein  
**Product:** Purified Recombinant Protein  
**Lot Number:**  
**Manufactured:**

**Catalog No.:** DP0012-P2  
**Concentration :**  
**EXP date:**

### Antibody Purification

**Starting Material:** cell culture supernatant  
**Buffer/Components:** 20mM PB, 150mM NaCl, 3mM EDTA, 1/1000 P300  
**Storage Conditions:** -20 °C, Avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Protein Concentration	BCA methods		
Purity	SDS-Page (reduced)	G3302	

Prepared By: \_\_\_\_\_  
Quality Control-Date

Reviewed By: \_\_\_\_\_  
Quality Assurance-Date