

Product Name: HIV-2 gp36 protein
Product: Purified Recombinant Protein
Lot Number:
Manufactured:

Catalog No.: DP0012-P4
Concentration :
EXP date:

Antibody Purification

Starting Material: cell culture supernatant
Buffer/Components: 20mM PB, 150mM NaCl, 3mM EDTA, 1/1000 P300
Storage Conditions: -20°C, Avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Protein Concentration	BCA methods		
Purity	SDS-Page (reduced)	G3302	

Prepared By: _____
Quality Control-Date

Reviewed By: _____
Quality Assurance-Date