



## Certificate of Analysis

**Product Name:** RSV B Protein  
**Product:** Purified Recombinant Protein  
**Lot Number:**  
**Manufactured:**  
**EXP date:**

**Catalog No.:** DP0037-P2  
**Concentration :**

### Antibody Purification

**Starting Material:** cell culture supernatant  
**Buffer/Components:** 20mM PB, 150mM NaCl, 3mM EDTA, 1/1000 P300  
**Storage Conditions:** -20°C, Avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Protein Concentration	BCA methods		
Purity	SDS-Page (reduced)	G3302	

**Prepared By:**

**Reviewed By:**

No.168 Mashang Meiliang Road | Wuxi , Jiangsu province, China | 0510-81830666

[www.origene.com](http://www.origene.com)

[custsupport@origene.com](mailto:custsupport@origene.com)