

**Product Name:** PLA2R1  
**Product:** Purified Recombinant Protein  
**Lot Number:**  
**Manufactured:**

**Catalog No.:** DP0049-P1  
**Concentration :**  
**EXP date:**

### Protein Purification

**Starting Material:** Culture supernatant  
**Buffer/Components:** PBS, 10% glycerol  
**Storage Conditions:** -80°C, avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Protein Concentration	Bradford protein assay method		
Purity	SDS-Page (reduced)	G3302	

Certified by: \_\_\_\_\_

Date: \_\_\_\_\_