

Product Name: APOH
Product: Purified Recombinant Protein
Lot Number:
Manufactured:

Catalog No.: DP0051-P1
Concentration :
EXP date:

Protein Purification

Starting Material: Culture supernatant
Buffer/Components: PBS, 10% glycerol
Storage Conditions: -80°C, avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Protein Concentration	Bradford protein assay method		
Purity	SDS-Page (reduced)	G3302	

Certified by: _____ **Date:** _____