

Product Name: PCNA
Product: Purified Recombinant Protein
Lot Number:
Manufactured:

Catalog No.: DP0052-P1
Concentration :
EXP date:

Protein Purification

Starting Material: Lysate
Buffer/Components: 20mM Tris-HCl,pH7.5,2mM EDTA,20% glycerol.
Storage Conditions: -20°C

| Test: | Method: | Document: | Results: |
|-----------------------|-----------------------|-----------|----------|
| Protein Concentration | Microplate BCA method | | |
| Purity | SDS-Page (reduced) | G3302 | |

Certified by: _____ **Date:** _____