

Product Name: MPO
Product: Purified Recombinant Protein
Lot Number:
Manufactured:

Catalog No.: DP0054-P1
Concentration :
EXP date:

Antibody Purification

Starting Material: Culture supernatant
Buffer/Components: PBS, 10%glycerol.
Storage Conditions: -80°C, avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Protein Concentration	BCA methods		
Purity	SDS-Page (reduced)	G3302	

Certified by: _____ **Date:** _____