



Certificate of Analysis

**Product Code:** S2276MA01-MA  
**Description:** Mouse Anti-IgE Purified Antibody  
**Clone:** CL206  
**Lot Number:**  
**Site of Manufacture:** Newark, DE

Antibody Purification

**Salt Cut:** Ammonium Sulfate Precipitation  
**Axial Flow Column:** Q-Sepharose  
**Presentation/Final Buffer:** 100mM Sodium Phosphate, pH 7.0  
**Storage Conditions:** 4°C (Range 2°C to 8°C)

Test:	Method:	Document:	Results:
Protein Concentration	A280	G3476	
Purity	HPLC	G3331	
Isoelectric Point	Isoelectric Focusing	G3305	
Isotype	Isostrip	G3346	

Animal Health Certification:

Mice used were healthy and specific pathogen free. The manufacturing process was performed with animals of U.S. origin.

**Prepared By:**

**Date Reviewed:**

**Reviewed By:**

**Date Reviewed:**