

Product Code: S2276MA02-MA
Description: Mouse Anti-IgE Purified Antibody
Clone: CL102.1
Lot Number:
Site of Manufacture: Newark, DE

Antibody Purification

Axial Flow Column: Protein A Sepharose, Fast Flow
Presentation/Final Buffer: 100mM Sodium Phosphate, pH 7.0
Storage Conditions: 4°C (Range 2°C to 8°C)

Test:	Method:	Document:	Results:
Protein Concentration	A280	G3476	
Purity	SDS-Page (reduced)	G3302	
Isoelectric Point	Isoelectric Focusing	G3305	
Isotype	ISOSTrip™	G3346	

Animal Health Certification:

Mice used were healthy and specific pathogen free. The manufacturing process was performed with animals of U.S. origin.

QA Reviewed By:

Date Reviewed: