



## Certificate of Analysis

**Product Name:** ACE2 mouse monoclonal antibody, clone OTI1G4  
**Product:** Purified Monoclonal Antibody  
**Lot Number:**  
**EXP Date:**

**Catalog No.:** TA804152  
**Concentration:**

### Antibody Purification

**Starting Material:** cell culture supernatant  
**Axial Flow Column:** Protein G Sepharose  
**Buffer/Components:** Freeze-dried in PBS, BSA/glycerol-free, reconstitute with PBS (pH7.3)  
**Storage Conditions:** -20°C, Avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOSTrip™	G3346	
Application	ELISA	Datasheet	

**Prepared By:**

**Reviewed By:**

9620 Medical Center Dr. | Rockville, MD20850 USA | +1 (888) 267-4436 | +1 (301) 340-3188

[www.origene.com](http://www.origene.com) [custsupport@origene.com](mailto:custsupport@origene.com)