

Product Name: SARS-CoV-2 N protein
Product: Purified Recombinant Protein
Lot Number:
Manufactured:

Catalog No.: TP750179
Concentration :
EXP date:

Antibody Purification

Starting Material: Lysate
Buffer/Components: PBS
Storage Conditions: -80°C

Test:	Method:	Document:	Results:
Protein Concentration	Bradford protein assay		
Purity	SDS-Page (reduced)	G3302	

Prepared By: _____
Quality Control-Date

Reviewed By: _____
Quality Assurance-Date