

Certificate of Analysis

Product Information

Product name: PLA2R1 protein

Catalog No: TP790153

Lot Number:

Quantity:

Concentration:

Buffer Formulation: PBS.

Storage Conditions: Store at - 80 °C. Avoid repeated freeze-thaw cycles. Stable for at least 3 years from receipt of products under proper storage and handling conditions.

Quality Control

Test Items	Specifications	Results:
Protein Concentration	Microplate BCA Protein Assay	Conformed
Purity(SDS-PAGE)	> __%	Conformed

Certified by: _____

Date: _____