

Product Name: UCH-L1 mouse monoclonal antibody, clone OTI8E3
Product: Purified Monoclonal Antibody **Catalog No.:** DM0035-B1
Lot Number:
EXP Date:

Antibody Purification

Starting Material: cell culture supernatant
Axial Flow Column: Protein G Sepharose
Buffer/Components: Freeze-dried in PBS, BSA/glycerol-free, reconstitute with PBS (pH7.3)
Storage Conditions: -20°C, Avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOSTrip™	G3346	
Application		Datasheet	

Prepared By: _____
Quality Control-Date

Reviewed By: _____
Quality Assurance-Date