

**Product Name:** UCH-L1 Rabbit Monoclonal Antibody, clone OTIR2F11  
**Product:** Purified Monoclonal Antibody  
**Lot Number:** **Catalog No.:** DM0035-S4  
**PD:** **Concentration:**  
**EXP:**

### Antibody Purification

**Starting Material:** cell culture supernatant  
**Axial Flow Column:** Protein G Sepharose  
**Buffer/Components:** Freeze-dried in PBS, BSA/Glycerol-free, reconstitute with PSB (pH7.3)  
**Storage Conditions:** -20°C, Avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Purity	SDS-Page (reduced)	G3302	
Application		Datasheet	

Prepared By: \_\_\_\_\_  
Quality Control-Date

Reviewed By: \_\_\_\_\_  
Quality Assurance-Date