

**Description:** Calibrator AHA Protein  
**Part number:** DP0005-P3  
**Lot Number:**

**Date of Manufacture:**

**Expiration Date:**

**Form:** Yellow powder

**Storage Conditions:** 2-8°C, Avoid repeated freeze-thaws.

**Method of application:** Recommended to use 1ml sterile water dissolve sample 6 hours before use.

Test:	Specification	result:	Pass/Fall:
Appearance	Yellow powder	Yellow powder	pass
Final pH*	6.0-8.0	7.0	pass
OD280*	Report		pass
Sodium Azida Concentration*	</=0.1% by spectrophotometric method	0.1%	pass
Protein Concentration*	AAT>600mg/dL		pass
	HPT>520mg/dL		pass
	AGP>400mg/dL		pass

\*Data from the liquid state.

**Prepared By:** \_\_\_\_\_  
 Quality Control-Date

**Reviewed By:** \_\_\_\_\_  
 Quality Assurance-Date