

Description: Calibrator AHA Protein
Part number: DP0005-P4
Lot Number:

Date of Manufacture:

Expiration Date:

Form: Yellow powder

Storage Conditions: 2-8°C, Avoid repeated freeze-thaws.

Method of application: Recommended to use 0.5ml sterile water dissolve sample 4 hours before use.

Test:	Specification	result:	Pass/Fall:
Appearance	Yellow powder	Yellow powder	pass
Final pH*	6.0-8.0	7.0	pass
OD280*	Report		pass
Sodium Azida Concentration*	</=0.1% by spectrophotometric method	0.1%	pass
Protein Concentration*	AAT>600mg/dL		pass
	HPT>520mg/dL		pass
	AGP>400mg/dL		pass

*Data from the liquid state.

Prepared By: _____
 Quality Control-Date

Reviewed By: _____
 Quality Assurance-Date