



Description: Anti-human Haptoglobin
Product Code: S1011G000-S5
Lot Number:
Stage/Grade: TIA Grade
Manufacture Date:
Source: Goat
Antigen Source: Normal Human Plasma
Expiration Date:
Storage: 2°C to 8°C

Test	Specification	Result	Pass/Fail
Appearance	Slightly yellow to amber color, essentially free from particulates		
Specificity	Monospecific precipitin reaction vs Normal human plasma and 2x concentrated human serum		
	Arc has similar electrophoretic mobility as Master Lot		
Final pH	pH 6.0-8.0		
Total Protein by A280	Report only		

Processing: Delipidated, defibrinated and filtered 0.2µm
Buffer: 20 mM Sodium Phosphate, 150 mM Sodium Chloride, 0.05%NaN₃
Country of Origin: USA
Hazard Testing and Warnings

Source animal health: Certified free of rinderpest, hoof and mouth disease, contagious bovine pleuropneumonia, and test negative for blue tongue.

Sodium Azide: Product contains sodium azide. Sodium azide may react with lead copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts" in the Manual guide-safety Management NO. CDC-22 issued by the Center for Disease and Control, Atlanta, GA 1976.

Intended Use: For Manufacturing, Processing and Relabeling.

Prepared by:

Reviewed by: