



Certificate of Analysis

Description: Anti-human C-Reactive Protein
Product Code: S3911G000-S5
Lot Number:
Stage/Grade: Serum/TIA
Expiration Date:
Source: Goat
Antigen Source: Normal Human Plasma/Serum
Date of Manufacture:
Storage: 2°C to 8°C

Test	Specification	Result	Pass/Fail
Appearance	Slightly yellow to amber color, essentially free from particulates		
Specificity	No precipitin arc vs. Normal human plasma and 2x Concentrated normal human serum		
	Monospecific precipitin arc vs. CRP material		
Total Protein	Report only		
Final pH	6.0 – 8.0		

Processing: Delipidated, defibrinated and filtered 0.2µm
Buffer: 0.02M Sodium Phosphate, 0.15M Sodium Chloride, 0.05% NaN₃
Country of Origin: USA
Hazard Testing and Warnings

Source animal health: Certified free of rinderpest, hoof and mouth disease, contagious bovine pleuropneumonia, and test negative for blue tongue.

Sodium Azide: Product contains sodium azide. Sodium azide may react with lead copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts" in the Manual guide-safety Management NO. CDC-22 issued by the Center for Disease and Control, Atlanta, GA 1976.

Intended Use: For Manufacturing, Processing and Relabeling.

Prepared by:

Reviewed by: