

**Description: Anti-Human C-Reactive Protein IgG Fraction (CRP-IgG)**  
**Product Code: S3911GND0-D0**  
**Lot Number:**  
**Stage/Grade: IgG Fraction**  
**Manufacture Date:**  
**Source: Goat**  
**Antigen Source: Normal Human Plasma/Serum**  
**Expiration Date:**  
**Storage: 2°C to 8°C**

Test	Specification	Result	Pass/Fail
Appearance	Essentially colorless and free of particulates		
Specificity	No precipitin arc vs. Normal human plasma and 3x Concentrated normal human serum		
	Monospecific precipitin arc vs. CRP material		
Final pH	pH 6.0-8.0		
Total Protein by A280	Report only		

**Processing:** Delipidated, defibrinated, Ammonium Sulfate precipitated and filtered 0.2µm  
**Buffer:** 20 mM Sodium Phosphate, 150 mM Sodium Chloride, 0.05%Na<sub>3</sub>  
**Country of Origin:** USA

#### **Hazard Testing and Warnings**

**Source animal health:** Certified free of rinderpest, hoof and mouth disease, contagious bovine pleuropneumonia, and test negative for blue tongue.

**Sodium Azide:** Product contains sodium azide. Sodium azide may react with lead copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts" in the Manual guide-safety Management NO. CDC-22 Issued by the Center for Disease and Control, Atlanta, GA 1976.

**Intended Use:** For Further Manufacturing or Research Use Only. Product meets European Union requirements for treated technical blood products. Not for Human or Animal Consumption.

Prepared by/Date: