

Certificate of Analysis

Production Information

Product name: S100B
Catalog No.: TP300277
Lot Number:
Quantity:
Concentration:
Buffer Formulation: 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol.
Tags: C-Myc/DDK*
Storage Conditions: Store at -80°C. Avoid repeated freeze-thaw cycles. Stable for at least 3 years from receipt of products under proper storage and handling conditions.

** DDK is the same as FLAG from Sigma*

Quality Control

Test Items	Specifications	Results
Protein Concentration	Microplate Bradford Protein Assay	
Purity(SDS-PAGE)	>80%	

Certified by: _____ Date: _____