

Certificate of Analysis

Product Information

Product name: GFAP

Catalog No: TP304548

Lot Number:

Quantity:

Concentration:

Buffer Formulation: 25mM Tris-HCl, pH 7.3, 100mM glycine, and 10% glycerol.

Tags: C-terminal Myc/DDK* tag

Storage Conditions: Store at -80°C . Avoid repeated freeze-thaw cycles. Stable for at least 1 year from receipt of products under proper storage and handling conditions.

Quality Control

Test Items	Specifications	Results:
Protein Concentration	Microplate Bradford Protein Assay	
Purity(SDS-PAGE)	>80%	

Certified by: _____

Date: _____