

Description: Anti-Human Apolipoprotein A1 (Apo A1)
Product Code: S1141G000-S5
Lot Number:
Stage/Grade: TIA
Manufacture Date:
Source: Goat
Antigen Source: Normal Human Plasma/Serum
Expiration Date:
Storage: 2°C to 8°C

Test	Specification	Result	Pass/Fail
Appearance	Slightly yellow to amber color, essentially free from particulates		Pass
Specificity	Monospecific precipitin reaction vs Normal human plasma and 2X concentrated human serum		Pass
	Arc has similar electrophoretic mobility as Reference Lot		Pass
Final pH	pH 6.0-8.0		Pass
Total Protein by A280	Report only		N/A

Processing: Delipidated, defibrinated and filtered 0.2µm
Buffer: 20 mM Sodium Phosphate, 150 mM Sodium Chloride, 0.05%NaN₃
Country of Origin: USA and/or meets EU Requirements for treated technical blood products

Hazard Testing and Warnings

Source animal health: Certified free of rinderpest, hoof and mouth disease, contagious bovine pleuropneumonia, and test negative for blue tongue.

Sodium Azide: Product contains sodium azide. Sodium azide may react with lead copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts" in the Manual guide-safety Management NO. CDC-22 issued by the Center for Disease and Control, Atlanta, GA 1976.

Intended Use: For Further Manufacturing or Research Use Only. Product meets European Union requirements for treated technical blood products. Not for Human or Animal Consumption.

Prepared by/Date:

Reviewed by/Date: