

Description: Anti-human B2M-IgG
Part Number: S3801GND0-D0
Lot Number:
Stage/Grade: IgG Fraction
Manufacture Date:
Source: Goat
Antigen Source: Human Urine
Expiration Date:
Storage: 2°C to 8°C

Test	Specification	Result	Pass/Fail
Appearance	Visibly free of particulates		
Specificity	Reacts with Human B2M		
Protein by A ₂₈₀	40 – 60 mg/mL		

Processing: Filtered 0.2µm
Buffer: 20mM Sodium Phosphate, 150 mM Sodium Chloride, 0.05% NaN₃
Country of Origin: USA

Hazard Testing and Warnings

Sodium Azide: Product contains sodium azide. Sodium azide may react with lead copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. For further information, refer to "Decontamination of Laboratory Sink Drains to Remove Azide Salts" in the Manual guide-safety Management NO. CDC-22 issued by the Center for Disease and Control, Atlanta, GA 1976.

Intended Use: For Further Manufacturing or Research Use Only. Product meets European Union requirements for treated technical blood products. Not for Human or Animal Consumption.

Prepared by/Date: