

**Product Name:** NP(Influenza A) mouse monoclonal antibody , clone 27A14  
**Product:** Purified Monoclonal Antibody **Catalog No.:** D678  
**Lot Number:** **Concentration:**  
**Production Date:**  
**Exp Date:**

### Antibody Purification

**Starting Material:** cell culture supernatant  
**Axial Flow Column:** Protein G Sepharose  
**Buffer/Components:** PBS (pH7.3), 0.1%NaN<sub>3</sub>  
**Storage Conditions:** Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	A280	G3476	
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOSTrip™	G3346	
Application	ELISA	Datasheet	

Prepared By: \_\_\_\_\_  
Quality Control-Date

Reviewed By: \_\_\_\_\_  
Quality Assurance-Date