

Product Name: NP(Influenza B) mouse monoclonal antibody , clone 31D24
Product: Purified Monoclonal Antibody **Catalog No.:** D688
Lot Number: **Concentration:**
Production Date:
Exp Date:

Antibody Purification

Starting Material: cell culture supernatant
Axial Flow Column: Protein G Sepharose
Buffer/Components: PBS (pH7.3), 0.1%NaN3
Storage Conditions: Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	A280	G3476	
Purity	SDS-Page (reduced)	G3302	> 95%
Isotype	ISOSTrip™	G3346	Mouse IgM
Application	ELISA	Datasheet	PASS

Prepared By: _____
Quality Control-Date

Reviewed By: _____
Quality Assurance-Date