

**Product Name:** sFlt1 mouse monoclonal antibody , clone 38C19  
**Product:** Purified Monoclonal Antibody **Catalog No.:** DM0058-B4  
**Lot Number:** **Concentration:**  
**Production Date:**  
**Exp Date:**

### Antibody Purification

**Starting Material:** cell culture supernatant  
**Axial Flow Column:** Protein G Sepharose  
**Buffer/Components:** PBS (pH7.3), 0.1%NaN3  
**Storage Conditions:** Store at 2-8°C ,Do Not freeze.

| Test:                 | Method:            | Document: | Results: |
|-----------------------|--------------------|-----------|----------|
| Protein Concentration | A280               | G3476     |          |
| Purity                | SDS-Page (reduced) | G3302     |          |
| Isotype               | ISOSTrip™          | G3346     |          |
| Application           | ELISA              | Datasheet |          |

Prepared By: \_\_\_\_\_  
Quality Control-Date

Reviewed By: \_\_\_\_\_  
Quality Assurance-Date