

Product Name: IgG rabbit monoclonal antibody , clone 200C6
Product: Purified Monoclonal Antibody **Catalog No.:** DM0078-B1
Lot Number: **Concentration:**
Production Date:
Exp Date:

Antibody Purification

Starting Material: cell culture supernatant
Axial Flow Column: Protein A Sepharose
Buffer/Components: PBS (pH7.3), 0.1%NaN3
Storage Conditions: Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	A280	G3476	
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOSTrip™	G3346	
Application	ELISA	Datasheet	

Prepared By: _____
Quality Control-Date

Reviewed By: _____
Quality Assurance-Date