

**Product Name:** ADV mouse monoclonal antibody , clone 27B15  
**Product:** Purified Monoclonal Antibody **Catalog No.:** DM0082-B2  
**Lot Number:** **Concentration:**  
**Production Date:**  
**Exp Date:**

### Antibody Purification

**Starting Material:** cell culture supernatant  
**Axial Flow Column:** Protein G Sepharose  
**Buffer/Components:** PBS (pH7.3), 0.1%NaN3  
**Storage Conditions:** Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	A280	G3476	
Purity	SDS-Page (reduced)	G3302	> 95%
Isotype	ISOSTrip™	G3346	Mouse IgG1
Application	ELISA	Datasheet	PASS

Prepared By: \_\_\_\_\_  
Quality Control-Date

Reviewed By: \_\_\_\_\_  
Quality Assurance-Date