

Product Name: NP FluA
Product: Purified Recombinant Protein
Lot Number:
Production Date:

Catalog No.: DP0013-P2
Concentration:
Expiring Date:

Protein Purification

Starting Material: Cell Lysate
Buffer/Components: 50 mM Tris-HCl, pH 8.0, 500 mM NaCl
Storage Conditions: Store at -80°C. Avoid repeated freeze-thaw cycles.

Test:	Method:	Document:	Results:
Protein Concentration	Microplate BCA method	G3476	
Purity	SDS-Page (reduced)	G3302	

Prepared By: _____
Quality Control-Date

Reviewed By: _____
Quality Assurance-Date