

**Product Name:** NP FluB  
**Product:** Purified Recombinant Protein  
**Lot Number:**  
**Production Date:**

**Catalog No.:** DP0014-P2  
**Concentration:**  
**Expiring Date:**

### Protein Purification

**Starting Material:** Cell Lysate  
**Buffer/Components:** 50 mM Tris-HCl, pH 8.0, 500 mM NaCl  
**Storage Conditions:** Store at -80°C. Avoid repeated freeze-thaw cycles.

Test:	Method:	Document:	Results:
Protein Concentration	Microplate BCA method	G3476	
Purity	SDS-Page (reduced)	G3302	

Prepared By: \_\_\_\_\_

Quality Control-Date

Reviewed By: \_\_\_\_\_

Quality Assurance-Date