

Product Name: PLA2R1**Product:** Purified Recombinant Protein**Lot Number:****Production Date:****Catalog No.:** DP0049-P1**Concentration:****Expiring Date:****Protein Purification****Starting Material:** cell culture supernatant**Buffer/Components:** PBS, 10% glycerol**Storage Conditions:** Store at -80°C. Avoid repeated freeze-thaw cycles.

Test:	Method:	Document:	Results:
Protein Concentration	Microplate BCA method	G3476	
Purity	SDS-Page (reduced)	G3302	

Prepared By: _____

Quality Control-Date

Reviewed By: _____

Quality Assurance-Date