

Product Name: PCNA-BIO
Product: Purified Recombinant Protein
Lot Number:
Production Date:

Catalog No.: DP0052-P2
Concentration:
Expiring Date:

Protein Purification

Starting Material: Lysate
Buffer/Components: 20mM Tris-HCl, pH7.5, 5mM EDTA, 20% glycerol
Storage Conditions: -80°C

Test:	Method:	Document:	Results:
Protein Concentration	Microplate BCA method		
Purity	SDS-Page (reduced)	G3302	

Prepared By: _____
Quality Control-Date

Reviewed By: _____
Quality Assurance-Date