

**Product Name:** MDA5  
**Product:** Purified Recombinant Protein  
**Lot Number:**  
**Production Date:**

**Catalog No.:** DP0053-P1  
**Concentration:**  
**Expiring Date:**

**Protein Purification**

**Starting Material:** Cell Lysate  
**Buffer/Components:** 25mM Tris-HCl, 100mM glycine, and 10% glycerol, pH 7.3  
**Storage Conditions:** -80°C

Test:	Method:	Document:	Results:
Protein Concentration	Bradford protein assay method		
Purity	SDS-Page (reduced)	G3302	

**Prepared By:** \_\_\_\_\_  
Quality Control-Date

**Reviewed By:** \_\_\_\_\_  
Quality Assurance-Date