

**Product Name:** MDA5-BIO  
**Product:** Purified Recombinant Protein  
**Lot Number:**  
**Production Date:**

**Catalog No.:** DP0053-P2  
**Concentration:**  
**Expiring Date:**

### Protein Purification

**Starting Material:** Cell Lysate  
**Buffer/Components:** 25mM Tris-HCl, 100mM glycine, and 10% glycerol, pH 8.0  
**Storage Conditions:** Store at -80°C. Avoid repeated freeze-thaw cycles.

Test:	Method:	Document:	Results:
Protein Concentration	Bradford protein assay method	G3476	
Purity	SDS-Page (reduced)	G3302	

**Prepared By:** \_\_\_\_\_  
Quality Control-Date

**Reviewed By:** \_\_\_\_\_  
Quality Assurance-Date