

**Product Name:** SEPSECS  
**Product:** Purified Recombinant Protein  
**Lot Number:**  
**Production Date:**

**Catalog No.:** DP0056-P1  
**Concentration:**  
**Expiring Date:**

**Protein Purification**

**Starting Material:** Cell Lysate  
**Buffer/Components:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol  
**Storage Conditions:** Store at -80°C. Avoid repeated freeze-thaw cycles.

<b>Test:</b>	<b>Method:</b>	<b>Document:</b>	<b>Results:</b>
Protein Concentration	Bradford protein assay method		
Purity	SDS-Page (reduced)	G3302	

**Prepared By:** \_\_\_\_\_  
Quality Control-Date

**Reviewed By:** \_\_\_\_\_  
Quality Assurance-Date