

**Product Name:** DSG1  
**Product:** Purified Recombinant Protein  
**Lot Number:**  
**Manufactured:**

**Catalog No.:** DP0057-P1  
**Concentration :**  
**EXP date:**

### Protein Purification

**Starting Material:** Cell Lysate  
**Buffer/Components:** 25mM Tris-HCl, 100mM glycine, 250mM NaCl, and 10% glycerol, pH 7.3.  
**Storage Conditions:** -80°C

Test:	Method:	Document:	Results:
Protein Concentration	Bradford protein assay method		
Purity	SDS-Page (reduced)	G3302	

**Certified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_