

**Product Name:** GAGE7  
**Product:** Purified Recombinant Protein  
**Lot Number:**  
**Production Date:**

**Catalog No.:** DP0059-P1  
**Concentration:**  
**Expiring Date:**

**Protein Purification**

**Starting Material:** Cell Lysate  
**Buffer/Components:** 25 mM Tris, 500 mM NaCl, 0.4 mg/ml Triple Flag, pH 7.3  
**Storage Conditions:** Store at -80°C. Avoid repeated freeze-thaw cycles.

<b>Test:</b>	<b>Method:</b>	<b>Document:</b>	<b>Results:</b>
Protein Concentration	Bradford protein assay method	G3476	
Purity	SDS-Page (reduced)	G3302	

**Prepared By:** \_\_\_\_\_  
Quality Control-Date

**Reviewed By:** \_\_\_\_\_  
Quality Assurance-Date